

Intro to Neuroradiology MRI

Luke Wojdyla, DO




Table of Contents

Night Float Protocolling Cheat Sheet

Brain/Head:

- "Brain WO"
 - Most patients with the indication of "stroke", "TIA", "AMS", "focal neurologic deficits"
- "Acute Stroke"
 - Only when the "neuro stroke" team orders it (most "stroke" indications should instead be "Brain WO" above). Abbreviated protocol for patients eligible for intervention. Ordering team coordinate with the MRI techs, potentially pulling other patients off the scanner. No discussion b/t stroke team and radiologist needed, but rads can call MRI to give them heads up
- "Brain Tumor"
 - Primary CNS neoplasms work ups, such as glioblastoma, medulloblastoma (not meningioma!)
- "Brain mets"
 - Intracranial metastases, infection
- Brain mass status post resection
 - For brain masses in the immediate post-operative period (24-72 hours), a full "Brain Tumor" protocol is not necessary. There is no prebuilt "epic" protocol button for this, but you can choose "Brain Mets" sequence and add any additional sequences requested by ordering (eg, T2 CUBE)
- "Ant Skull Base + Sella"
 - Choose for all "MRI Brain + Sella" orders, all pre-op or immediate post-op studies for sella or anterior skull base masses.
 - Note: "MRI Sella" protocol is a distinct protocol appropriate for Brain and Sella orders to evaluate for pituitary microadenoma or cases that are NOT for preoperative evaluation.
 - Pituitary apoplexy should be MRI Sella
- "MS Brain" +/- contrast
 - Suspected demyelination (MS, NMO, MOGAD, ADEM)
 - MRI brain, orbits, cervical spine, and thoracic spine may all be appropriate depending on scenario. If concerned for an initial or acute (exacerbation), contrast may be appropriate. Can generally just defer these to whatever neurology wants at the given time.

Spine

- Contrast not needed for trauma or degenerative myelopathy. Contrast is useful for infection, sometimes when evaluating the post-operative spine (complications), neoplasm, and sometimes demyelination (MS).
- Non-contrast spines are often auto-protocolled, but techs may sometimes still ask you for guidance. In cervical spine there are 2 trauma protocols, one called a "CCJ" (craniocervical junction) and the other "trauma". "Cervical CCJ" offers a dedicated look at the occiput-C1-C2 complex, specifically to assess for ligamentous injuries, but takes longer to do. "Cervical Trauma" can be used if injury on CT doesn't include CCJ,
- "Spine Survey" (a modified combined Cervical, Thoracic, and Lumbar spine WO/W) should be used for metastases, non-localized infection, and acute myelopathy/cord compression. Not infrequently, ortho/neurosurg wants a dedicated look at one section (eg, dedicated lumbar spine) AND screen the rest of the spine. In this case, recommend speaking with technologists, can adjust field of view on spine survey so that the entire Cervical and Thoracic spine are included in the upper portion, and the lower portion is just the lumbar spine with added pre-contrast axial images (for lumbar spine, ax T2 and ax T1). Or can just do dedicated Cervical, Thoracic, and Lumbar spine WO/W, but many patients will struggle tolerating being in the scanner that long

MRA:

- For Brain, non-contrast are sufficient for intracranial aneurysm screening and LVO detection (eg, time of flight and phase-contrast techniques)
- For neck, specifically add "T1 fat-saturated black blood" sequence (a non-contrast technique) to assess for dissection
- Contrast useful for more accurate vessel lumen assessment (especially in the neck) for high-grade carotid stenosis, dissections, and aortic arch pathology, and can also be useful to characterize treated/coiled aneurysms to look for residual filling
- For aneurysms after surgical clipping, CTA is superior.
- For aneurysms after endovascular coiling, MRA is superior.

MRV

- For venous sinus thrombosis, contrast-enhanced is **strongly** preferred (push for this). Exception is pregnant patients, can only do non-contrast.

Table of Contents

Techniques/Abbreviations

<u>1.5 vs 3 T</u>	<u>FLASH: Fast Low Angle Shot (FLASH, Siemens)</u>
<u>ASL: Arterial Spin Labeling</u>	<u>InHance 3D Velocity: 3D phase contrast MR angiography</u>
<u>Black Blood T1 Fat-Saturated (dissection protocol)</u>	<u>LAVA: Liver Acquisition with Volume Acceleration</u>
<u>BRAVO: Brain Volume Imaging</u>	<u>MARS: Metal Artifact Reduction Sequences</u>
<u>Brainlab: Surgical navigation system</u>	<u>MSDE: Motion-Sensitized Driven Equilibrium</u>
<u>Chemical Shift Fat Suppression</u>	<u>MERGE: Multi-Echo Recombined Gradient Echo</u>
<u>CISS (Constructive Interference in Steady State)</u>	<u>MP-RAGE: Magnetization-Prepared Rapid Acquisition Gradient Echo</u>
<u>Contrast Enhanced Techniques</u>	<u>MUSE: Multiplexed Sensitivity-Encoding</u>
<u>CUBE: Volumetric Isotropic Fast Spin Echo</u>	<u>PROPELLER (PROP): Periodically Rotated Overlapping Parallel Lines with Enhanced Reconstruction</u>
<u>DSC: Dynamic susceptibility contrast</u>	<u>SPACE: Sampling Perfection with Application-optimized Contrasts using different flip angle Evolutions</u>
<u>DWI: Diffusion-Weighted Imaging</u>	<u>STIR: Short TI Inversion Recovery</u>
<u>FIESTA: Fast Imaging Employing Steady-State Acquisition</u>	<u>SWAN: Susceptibility-Weighted Angiography</u>
<u>FLAIR: Fluid Attenuated Inversion Recovery</u>	<u>TIRM: Turbo Inversion Recovery Magnitude</u>
<u>FLEX: Fast Leveraged Excitation</u>	<u>TOF: Time of Flight</u>
<u>FrFSE: Fast Recovery Fast Spin Echo</u>	<u>TRICKS: Time-Resolved Imaging of Contrast Kinetics</u>
<u>FSE: Fast Spin Echo</u>	<u>VIBE: Volumetric Interpolated Breath-hold Examination</u>
<u>FSPGR: Fast Spoiled Gradient-Recalled</u>	

UCSD Magnets

- ACTRI
 - 3T GE Signa
- CTIPM1 MR00W4
 - 3T GE Signa
- EC
 - 1.5T Siemens Symphony
- El Centro Regional Medical Center
 - 1.5T GE Artist
 - 1.5T GE HDxt
- HC MR3
 - 1.5T GE Signa (HDXT?)
- HC MR4
 - 3T GE Discovery
- HC HCOPCANMR1
 - 3T Canon Galon
- McG OPP 1
 - 3T Siemens Magnetom Vida
- McG OPP 2
 - 3T Siemens Magnetom Vida
- KOP MR1
 - 3T GE Signa Premier
- KOPMR2
 - 3T GE Signa Premier
- La Jolla (URLJ)
 - 1.5 T GE Signa HDXT
- TH TMR
 - 1.5T GE Signa (HDXT?)
- TH Mobile
 - 1.5T Siemens Magnetom Aera

UCSD Neuroradiology MRI protocols

[Brain](#)

- [Acute Stroke](#)
- [Brain WO](#)
- [Brain +C](#)
- [Brain Mets](#)
- [Brain Tumor](#)
- [Brain epilepsy W/WO](#)
- [Ant Skull Base + Sella](#)
- [MS W/WO contrast](#)

[Vascular](#)

- [Head MRA](#)
- [Head MRV](#)
- [Neck MRA](#)
- [Spine MRA](#)

[Head and Neck](#)

- [Face](#)
- [Neck](#)
- [Sella/Pituitary](#)
- [Orbits](#)
- [IAC W/WO](#)
- [Trigeminal](#)
- [Glossopharyngeal](#)

[Spine](#)

- [C-spine WO](#)
- [C-spine CCJ](#)
- [C-spine WO/W](#)
- [C-spine MS WO/W](#)
- [T-spine WO](#)
- [T-spine Trauma](#)
- [T-spine MS WO/W](#)
- [T-spine WO/W](#)
- [L-spine WO](#)
- [L-spine trauma](#)
- [L-spine WO/W](#)
- [Spine survey](#)

Epic button pending

- Operative planning (brain)
- CSF leak (spine)
- AVM (brain)

Techniques

1.5T vs 3T MRI

3T



- Advantages
 - Higher signal-to-noise ratio (SNR)
 - Better spatial resolution or faster scan times
 - Improved lesion conspicuity (e.g., small MS plaques, hippocampal detail, subtle cortical dysplasia)
 - Better advanced imaging: ASL, DTI, fMRI, spectroscopy, and high-resolution vessel wall imaging
- Disadvantages
 - More susceptible to artifacts:
 - Susceptibility (e.g., near sinuses/skull base)
 - Chemical shift, dielectric shading
 - SAR (Specific Absorption Rate) limits may restrict sequence parameters (esp. in larger patients)
 - CSF flow/pulsation artifacts are more pronounced in spine imaging
 - Implant compatibility more restrictive

1.5T


- Advantages
 - Fewer artifacts in areas prone to susceptibility (e.g., skull base, posterior fossa)
 - Better fat suppression consistency (especially with chemical fat sat)
 - Less SAR concern → greater flexibility in sequence design
 - More implant-compatible
 - Spine imaging generally more robust in challenging regions (e.g., thoracic spine)
- Disadvantages
 - Lower SNR → may need longer scan times
 - Less sensitive for small lesions, especially in brain (e.g., cortical/subcortical MS, early infarcts)
 - Lower resolution for advanced techniques (e.g., DTI tractography, ASL perfusion, MR spectroscopy)

ASL


- ASL (Arterial Spin Labeling)
 - Non-contrast perfusion technique that uses RF pulses to magnetically label inflowing arterial blood as an endogenous tracer. Images are subtracted from a baseline control (no label) to quantify cerebral blood flow (CBF)
 - Used in stroke, tumor perfusion, and dAVF detection, can show hyperperfusion in AV shunting or hyperemia
 - 3D FSE readout with stack-of-spirals
 - Longer acquisition (4–6 min)




CCF



Chronic venous thrombosis



Hemangioblastoma




Jugular venous reflux

[Back to techniques](#)



Black blood T1 fat-saturated (MRA)

- Dissection protocol (primarily in the neck)
 - Most black blood T1 sequences are **fast spin echo (FSE/TSE)**–based
 - They use **double inversion recovery (DIR)** or **spatial presaturation bands** to suppress blood signal ("black blood")
 - Chemical shift fat saturation with a **frequency-selective saturation pulse** that targets the fat resonance
- Suppressing fat makes it easier to detect:
 - **Intraluminal thrombus**
 - **Intramural hematoma** in vessel walls (which appear T1 hyperintense)
 - **Wall thickening** or enhancement



Brainlab

- Surgical navigation and planning system that integrates with preoperative imaging (MRI, CT)
Used to define surgical targets, trajectories, and margins
- Supports intraoperative neuronavigation
- Relies on high-resolution 3D imaging (e.g., post-contrast MP-RAGE, FSPGR, or CUBE) with complete head coverage including the nose for registration
- Often used in brain tumor resections, epilepsy surgery, and functional neurosurgery
- Some institutions use StealthStation (Medtronic) as an alternative system



[Back to techniques](#)


Chemical Shift Fat Suppression

- Uses a short radiofrequency pulse tuned to the Larmor frequency of fat protons resulting in selective excitation followed by gradient spoiling, leaving only water protons to contribute to signal when the main excitation pulse is delivered shortly thereafter
- **Inhomogenous fat suppression** further from isocenter and around sources of susceptibility
- Compared to STIR, chemical shift fat suppression has higher SNR and is faster



CUBE (GE) or SPACE (Siemens)

- 3D volumetric TSE high spatial resolution and sub-millimeter isotropic voxels
- Contrast mechanisms differ from 2D TSE due to variable flip angle refocusing pulses, phase-encoding strategy
- Can be T1-, T2-, or FLAIR-weighted
- Can be combined with fat suppression, post-contrast
- SPACE: Sampling Perfection with Application-optimized Contrasts using different flip angle Evolutions
- Advantages
 - High spatial resolution (small structures like cranial nerves)
 - Operative planning
 - Enables high-quality multiplanar reformats
- Disadvantages
 - Lower SNR compared with 2D T2 TSE
 - Higher scan time compared to 2D sequences
 - Longer acquisition time (prone to motion artifact, can be accelerated with parallel imaging techniques)
 - CSF pulsation artifact



Whole brain 3D T2 CUBE


FIESTA

- Contrast is T2/T1-ratio-weighted, **not pure T2**
- Slightly higher **spatial resolution**
- More CSF pulsation artifact
- Faster acquisition



T2 CUBE

- Contrast is true T2 weighted
- Better T2 weighting and **tissue contrast**
- Less CSF pulsation artifact
- Slower acquisition


FIESTA




T2 CUBE




- Dynamic susceptibility contrast
 - A dynamic T2*-weighted imaging technique used to assess cerebral perfusion
 - A single bolus of gadolinium contrast is rapidly injected (4–5 mL/sec)
 - As the bolus passes through the cerebral vasculature, a T2*-weighted echo planar imaging (EPI) sequence (usually every 1–2 seconds) captures the signal drop caused by susceptibility effects of gadolinium as it passes through the brain's capillary network
 - Generates quantitative perfusion maps including CBF (cerebral blood flow), CBV (cerebral blood volume), and MTT (mean transit time)
 - Requires arterial input function (AIF) selection and post-processing for map generation
 - Commonly used in stroke, tumor, and vascular imaging to evaluate perfusion deficits or hypervascularity
 - Sensitive to susceptibility artifacts (given gradient echo EPI readout)




Perfusion
(source)



rCBV




rCBF



MTT

DWI (Diffusion-Weighted Imaging)

- Measures the random (Brownian) motion of water molecules in 3 orthogonal dimensions
 - Primarily used in stroke imaging, tumor characterization, and infection
 - Single-shot spin-echo echo-planar imaging (EPI) technique
 - Spin-echo core provides diffusion contrast
 - EPI readout is gradient-echo-based, giving inherent T2* susceptibility sensitivity (can detect hemorrhage)
 - Acquires all k-space lines from one RF excitation per slice (single TR) → very fast
 - Two sets of images from the same acquisition, differing in diffusion-gradient strength (b-value):
 - Higher b values represent greater diffusion weighting relative to T2 weighting
 - Degree of diffusion weighting is determined by the amplitude and duration of applied diffusion gradients
- b0
 - Diffusion gradients off
 - Essentially a T2-weighted image. Useful to troubleshoot when other sequences are motion-degraded or not acquired (e.g., missing SWI)
 - b1000 (AKA "trace," labeled "DWI" in PACS)
 - Same acquisition with strong diffusion gradients applied
 - ADC (apparent diffusion coefficient)
 - Calculated, not separately acquired; a quantitative measure of diffusion (removes T2 shine-through)
 - eADC (exponential ADC)
 - Further derived from ADC; qualitative, retains DWI-like appearance (restriction stays bright)



DWI (Diffusion-Weighted Imaging)

- Pros


- Extremely fast (single-shot) → motion-robust, seconds per acquisition
- Most sensitive sequence for acute infarct — positive within minutes
- Aids tumor grading, abscess vs necrotic tumor, cholesteatoma, epidermoid
- b_0 doubles as a T_2^* -sensitive backup to screen for blood products

- Limitations

- Low spatial resolution and inherent image blur
- Poor soft-tissue anatomic detail vs conventional sequences
- Geometric distortion from the EPI readout
- Sensitive to susceptibility (air–bone interfaces, metal, hemorrhage)

- Pitfalls

- T2 shine-through
 - Hyperintensity on DWI can represent true diffusion restriction, intrinsic T_2 -hyperintensity also contributes to DWI signal intensity (T2 shine-through) and should not be confused for diffusion restriction
 - True diffusion restriction will exhibit DWI hyperintensity and corresponding ADC hypointensity -always confirm with low ADC
- T2 blackout
 - Very low T_2 signal can mask true restriction on b_{1000}
- Skull base / posterior fossa dropout and distortion near air (temporal bone, sinuses)
- ADC is quantitative but absolute values vary with b-value and vendor/technique



[Back to techniques](#)

DWI - Variants


- **PROPELLER**

- Non-EPI readout (fast-spin echo), acquires multiple overlapping radial "blades" in k-space and averages them for artifact reduction.
- Results in both higher resolution and significantly reduced motion artifacts
- Decreased FOV offsets the increase in acquisition time
- Lower susceptibility artifacts than ss-EPI (less distortion in air-bone interfaces like the skull base)
- Higher SNR than standard DWI due to multiple overlapping k-space acquisitions.
- Longer scan times than MUSE


- **MUSE (Multiplexed Sensitivity-Encoding)**

- Multi-shot EPI sequence designed to improve resolution and reduce distortion compared to standard single-shot EP
- Acquires data over multiple shots rather than a single shot
- Higher resolution
- Reduced geometric distortion (better for regions affected by susceptibility artifacts, such as the skull base and orbits)
- Longer scan times than standard, more motion artifact


DWI MUSE



DWI PROP




DWI (standard SS -EPI)



Recurrent cholesteatoma
only visible on PROP

[Back to techniques](#)




Single Shot

Multi-Shot

FIESTA (GE) or CISS (Siemens)

- FIESTA: Fast Imaging Employing Steady-State Acquisition
- CISS (Constructive Interference in Steady State)
- Balanced Steady-state free precession (SSFP) sequence
 - Runs at very short TR/TE with a steady-state maintained by repeated RF pulses.
 - Contrast is T2/T1-ratio weighted, not true T2 — fluid and fat are bright because of their high T2/T1 ratio



- Advantages
 - Highest SNR efficiency of any sequence (SNR per $\sqrt{\text{time}}$), signal recycled every TR
 - Very fast, motion-robust
 - Crisp delineation of small structures bathed in CSF
 - 3D isotropic
- Disadvantages
 - Banding artifacts — dark off-resonance bands, worsening at high field and with imperfect shim.
 - Susceptibility-sensitive at air–tissue interfaces (skull base, sinuses).
 - Not true T2 — poor soft-tissue-to-soft-tissue contrast, weak for solid lesion characterization.
 - Fat is bright which can obscure pathology without fat suppression


FIESTA

- Contrast is T2/T1-ratio-weighted, **not pure T2**
- Slightly higher **spatial resolution**
- More CSF pulsation artifact
- Faster acquisition



T2 CUBE

- Contrast is true T2 weighted
- Better T2 weighting and **tissue contrast**
- Less CSF pulsation artifact
- Slower acquisition

FIESTA





T2 CUBE




FLAIR


- FLAIR (Fluid Attenuated Inversion Recovery)
 - 180 RF pulse suppresses free fluid signal (eg, CSF) at time T1 (inversion time)
 - Enhances visualization of pathology like white matter lesions
 - T1 FLAIR sequences used at most 3T scanners, some sequences on 1.5 T (e.g., on 1.5T, axial T1 incorporates FLAIR, but not sagittal)




Sag T1 FLAIR



Ax T2 FS FLAIR




Ax T1 FLAIR



3D T2 FS FLAIR C+

FLEX

- FLEX (Fast Leveraged EXcitation)
 - GE's implementation of the Dixon two-point dual-echo Dixon-based fat-water separation, producing in-phase, out-of-phase, water-only, and fat-only images
 - Can be added to various sequences
 - T1 or T2 FSE
 - GRE sequences ([LAVA](#), TRICKS)
 - Usually replaces STIR, post-contrast T1 sequences
 - Sequence nomenclature in PACS
 - InPhase = Standard T1 or T2 weighting
 - Water only image (fat suppressed, similar to STIR)
 - Fat only image (water suppressed, can be used for trouble-shooting)




Sag T2 FLEX (Water)



Sag T2 FLEX (InPhase)

FLEX

- While FLEX (Dixon-based fat-water separation) is excellent for fat suppression, it can occasionally cause **misregistration**, **fat-water swapping**, or **inhomogeneous suppression** in complex anatomy or with motion.
- The Fat-only sequences (fat is hyperintense) are often automatically sent to PACS for trouble-shooting purposes in case this occurs




FSE/TSE vs FrFSE/RST

- Fast Spin Echo (GE) or Turbo Spin Echo (Siemens)
 - "Default" for 2D T1 and T2 spin echo sequence, simply labeled as "Ax T1" or "Cor T2" (most 3T scanners)
 - Spin-echo sequence that acquires multiple echoes per TR, reducing scan time. After the last echo in the echo train, the spins must naturally recover back to equilibrium before the next TR. This leads to some signal loss between TRs, reducing overall image contrast for long-T2 tissues like CSF, edema, or joint fluid
- Fast Recovery Fast Spin Echo (GE) or Restore (Siemens)
 - FSE sequence with reduced signal loss due to improved recovery of longitudinal magnetization
 - A 90° fast recovery pulse is applied immediately after the last echo, forcing the remaining transverse magnetization back into the longitudinal direction. This shortens the T1 recovery time and allows the tissue to reach full magnetization faster before the next TR. As a result, fluid retains more signal and may appear brighter than in standard FSE.
 - Used in 2D T2 spin-echo sequences at 1.5T scanners

AX T2
09:23:28

FA: 142
TR: 6826 ms
TE: 113.7 ms
SE 0\512\352\0
NEX: 1
ARDL: 0.75\High

Series 6 IMA 20




pos/22.5 mm
thk/4 mm
sp/4.5 mm
FoV: 220 x 220 mm²
PB: 177.559

FSE

Ax T2 FrFSE
21:24:56

FA: 111
TR: 6771 ms
TE: 120.64 ms
SE 0\440\320\0
NEX: 1
ARDL: 0.75\High

Series 6 IMA 19




pos/-6.4 mm
thk/4 mm
sp/4.5 mm
FoV: 240 x 240 mm²
PB: 162.773

FrFSE

InHance 3D Velocity


- InHance 3D Velocity
 - A non-contrast, flow-sensitive MRA technique
 - Volumetric 3D phase contrast acquisition using velocity-encoding gradients to measure flowing blood
 - Generates magnitude and phase images
 - Can be reformatted into MRA or MRV projections
 - Best for: slower flowing structures, can evaluate venous structures on MRA



Degree of phase shift after gradient application indicates velocity

MARS (Metal artifact Reduction Sequences)

- Susceptibility mismatch between metal and tissue distorts B_0 → signal voids, bright "pile-up," geometric distortion (in-plane and through-plane), and failed fat suppression
- Basic strategies
 - Prefer 1.5T over 3T — artifact scales with field strength
 - Use TSE/FSE, avoid GRE — 180° refocusing pulses recover susceptibility dephasing
 - Prefer STIR over chemical-shift fat-sat (fails in the inhomogeneous field near metal; STIR's inversion recovery is field-robust)
 - Swap phase/frequency direction to reposition artifact away from the cord/region of interest
 - ↑ Receiver bandwidth
 - ↓ TE, ↑ echo train length, smaller voxels/thinner slices — reduces dephasing and through-plane blur
- Dedicated MAR sequences VAT (view-angle tilting)
 - Corrects in-plane distortion
- Tradeoffs / pitfalls
 - MAR sequences cost scan time, SNR, and/or spatial resolution
 - Over-correction can blur anatomy; residual distortion still possible with large hardware




T2


T2 MARS

MDSE Motion-Sensitized Driven Equilibrium



- Blood-suppression ("black-blood") nulls luminal blood so arterial wall can be assessed for atherosclerotic plaque, vasculitis, dissection, and aneurysms
- Mechanism
 - A preparatory module consisting of a 90° tip-down pulse $\rightarrow$ diffusion/motion-sensitizing gradients $\rightarrow$ 180° refocusing pulse $\rightarrow$ gradients $\rightarrow$ 90° tip-up pulse
 - Flowing blood accumulates phase dispersion from the motion-sensitizing gradients and is driven to zero signal (black blood), while stationary tissue (vessel wall) is unaffected and retains signal
- Pros
 - Blood suppression is independent of flow velocity/direction — works for slow, recirculating, or in-plane flow
 - Uniform suppression over large 3D volumes
- Limitations
 - SNR penalty — preparation module introduces T2 decay during the MSDE interval, reducing overall vessel wall SNR
 - Eddy currents from the motion-sensitizing gradients $\rightarrow$ incomplete restoration, artifacts
 - Longer scan time from the added prep
 - Over-strong gradients suppress the wall itself, not just blood
 - Incomplete CSF suppression can mimic or obscure outer-wall enhancement
 - Residual slow-flow signal can mimic wall thickening or plaque



3D MSDE





3D MSDE C+



3D MSDE C+

MERGE



- Multi-Echo Recombined Gradient Echo
 - A multi-echo GRE sequence
 - Combines echoes to improve SNR and contrast for evaluating vertebral endplates, disc margins, and subtle marrow pathology
 - T2* weighting enhances susceptibility effects makes it useful for detecting hemorrhage (will appear hypointense)



T2 hypointense structure (osteophyte vs disc) narrows **neuroforamina** or **spinal canal**, which is hyperintense on MERGE, confirming disc (not osteophyte)



Hypointense epidural hematoma



T1 post-contrast: spin echo (2D) versus MP-RAGE (3D)




Spin echo (2D)

Strengths

- High contrast to noise (CNR)
- Fast

Weaknesses

- Thick slices with gaps
- Pulsation artifact



[Back to techniques](#) T1 FSE C+



MP-RAGE (3D)

Strengths

- High signal to noise (SNR)
- High spatial resolution
- Surgical navigation
- Tiny lesions

Weaknesses


- Lower contrast per voxel
- More sensitive to motion
- Slow



MP-RAGE C+

MP-RAGE/BRAVO


- Magnetization-Prepared Rapid Acquisition Gradient Echo (MP-RAGE, Siemens) /Brain Volume Imaging (BRAVO, GE)
 - 3D T1-weighted GRE sequence with an 180° inversion-recovery prep pulse before the GRE readout train (increasing stronger longitudinal T1 contrast)
 - Gray-white matter differentiation superior to standard T1 FSE enhance gray-white matter contrast
 - Used in structural brain imaging and pre/post-contrast studies.
 - Better contrast resolution than FSPGR/BRAVO
 - IR prep can partially suppress T1 enhancement (correlate with FSE)



[Back to techniques](#) 2D T1 FSE C+

FSPGR/FLASH

- Fast Spoiled Gradient Echo (FSPGR, GE) / Fast Low Angle Shot (FLASH, Siemens)
- Similar to MP-RAGE, but no inversion pulse
 - T1 contrast comes from short TR + low flip angle alone
 - No IR results in decreased tissue contrast, but enhancement isn't suppressed




3D MP-RAGE C+

Some lesions are only seen on 2D T1 spin echo



2D T1 C+




2D T1 C+



3D MPRAGE C+

...some lesions are only seen on MP-RAGE



2D T1 C+




3D MP-RAGE C+


Must look at both sequences!

LAVA/VIBE


- Liver Acquisition with Volume Acceleration (LAVA, GE) and Volumetric Interpolated Breath-hold Examination (VIBE, Siemens)
- Mechanism
 - 3D gradient-echo T1-weighted sequence that uses the Dixon technique (see [FLEX](#)) to acquire in-phase and out-of-phase images in a single breath-hold
- Uses
 - Post-contrast sequences in spine, neck, body
 - Only T1-weighted, not applicable for T2 weighted images
- Pros:
 - Thin isotropic 3D
 - Fat-suppressed
 - Fast enough for breath-hold and multiphase dynamic post-contrast
- Limitations:
 - Interpolation blur
 - Lower SNR per voxel
 - Fat-sat failure
 - Susceptibility near air/metal/interfaces




Water: Ax LAVA T1 FS C+



InPhase: Ax LAVA T1 FS C+



Water: Cor LAVA T1 FS C+




Water: Cor LAVA T1 FS

Post-contrast Acquisition Timing


- T1 C+:
 - Peak T1-shortening effects of gadolinium occurs 5–10 minutes after injection
 - Imaging should be performed at peak enhancement
- T2-FLAIR C+:
 - Performed immediately after gadolinium administration while waiting for peak enhancement on T1+C sequences

Post-contrast Brain protocols


Brain C+ Default



Cor T1 FS +C




Ax T1 +C




T1 MP-RAGE +C


Brain Mets Leptomeningeal
(infection, mets)




Cor T1 FS +C



Ax T1 +C




T1 MP-RAGE +C




T2 FLAIR CUBE +C


Brain Tumor Primary* brain tumor
(glioma)



Ax T1 +C



T1 MP-RAGE +C




DSC Perfusion C+

[Back to techniques](#)

PROPELLER

- Periodically Rotated Overlapping Parallel Lines with Enhanced Reconstruction
 - Motion-resistant multi-shot fast spin-echo
 - BLADE for Siemens
 - Acquires k-space in rotating blades, reducing motion and pulsation artifacts
 - Used in both T1 and T2 FSE, DWI
 - Scan time longer than standard




Cor T2




Cor T2 PROP

STIR

- STIR (Short TI Inversion Recovery):
 - TIRM (Turbo Inversion Recovery Magnitude) on Siemens
 - Fat-suppressed T2-weighted inversion recovery sequence
 - Cannot be used post-contrast (nulls enhancement)
- STIR vs. T2 spin echo:
 - Fat Suppression Method:
 - STIR: Uses inversion recovery (short TI) to null fat signal based on T1 relaxation.
 - T2 Fat-Sat: Uses frequency-selective fat saturation pulses based on chemical shift difference between fat and water.
 - Field Inhomogeneity Sensitivity:
 - STIR: Insensitive to magnetic field inhomogeneity — performs reliably near metal, air-tissue interfaces, or implants.
 - T2 Fat-Sat: Susceptible to field inhomogeneity — fat suppression may fail near hardware or in areas with poor shimming (e.g., skull base, sinuses).
 - Signal-to-Noise Ratio:
 - STIR:
 - Lower SNR compared to spine echo
 - Contrast-to-noise for certain lesions may be higher
 - T2
 - Generally higher SNR
 - Speed
 - STIR takes longer, T2 is faster



At short to medium TI's, long T1 substances (such as CSF) are bright on magnitude IR. This is opposite to SE where long T1 values cause reduced signal intensity.



Additive T1 and T2 effects. MS plaque seen better on STIR image (right) than SE image (left)

On STIR, even if the overall image is grainier (lower SNR), the lesion is easier to detect because of the strong signal difference relative to background.

SWAN

- Susceptibility-Weighted Angiography
 - 3D multi-echo GRE sequence with magnitude and filtered phase combination
 - SWI on Siemens
 - Highlights susceptibility effects from blood products (paramagnetic), iron and calcium (diamagnetic)
 - Used for microbleeds, hemorrhage, thrombosis, calcification, and venous structures

Table 1: Overview of Susceptibility-Sensitive Sequences for MRI Vendors


Handedness	Sequence Name	Paramagnetic: Probable Microbleed	Diamagnetic: Probable Microcalcification
Left-handed		Hyperintense on phase images	Hypointense on phase images
Siemens	SWI		
Canon	FSBB		
Philips	3D GRE raw data		
Right-handed		Hypointense on phase images	Hyperintense on phase images
United Neusoft	SWI		
GE Healthcare	SWAN		
Philips	3D SWIp		

Note.— FSBB = flow sensitive black blood, GRE = gradient echo, SWAN = susceptibility-weighted angiography, SWI = susceptibility-weighted imaging, SWIp = SWI with phase enhancement, 3D = three-dimensional.


SWAN - Interpretation

- Step 1: Determine handedness of scanner
 - GE scanners are right-handed, so paramagnetic substances (blood) will appear hypointense on phase map, diamagnetic substances (calcium) will appear hyperintense
 - Use **internal cerebral veins** as an internal control to verify
- Step 2: If there is **both hyperintense and hypointense** signal present, use the **dipole effect** to troubleshoot
 - Dipole effect:
 - Predictable blooming artifact (magnetic field distortion pattern) that extends along the direction of B_0 (craniocaudal) caused by a localized susceptibility source
 - Diagram (right) is for right-handed scanners (GE), would be opposite for left-handed scanners (Siemens)


Calcification
(diamagnetic)




Blood
(paramagnetic)




Confirmed bleed (paramagnetic)




SWAN




Phase



Slice above



Slice below




Coronal

TOF

- Time of Flight
 - A non-contrast technique that visualizes flow-related enhancement from inflowing unsaturated blood
 - T1 weighted 2D or 3D gradient echo with short TR/TE and high flip angle (such as FSPGR)
 - Suffers from saturation in slow flow or in-plane flow
 - 3D MOTSA (Multiple Overlapping Thin Slab Acquisition) often used for neck


2D TOF

- Fine slices are stacked alongside each other
- Sensitive to slow flows (which do not remain in the slice for long and will therefore not be saturated)
- Poor spatial resolution along the axis of the slice stack




3D TOF

- Good spatial resolution in all 3 planes
- Better signal-to-noise ratio
- Longer acquisition




TOF

- Stationary tissue becomes saturated by repeated RF pulses
- Blood flow perpendicular to imaging plane is not saturated and has high signal ("flow-related enhancement")




Adapted from <https://doi.org/10.1148/rg.2021210141>



[Back to techniques](#)

TRICKS


- Time-Resolved Imaging of Contrast Kinetics
- 3D contrast-enhanced acquisition by sampling central k-space (AKA keyhole) more frequently than peripheral k-space
- Generates multiple time-resolved frames during contrast bolus passage Useful for evaluating vascular malformations, dAVFs, and high-flow lesions
- Can be applied to brain, spine, and peripheral vascular imaging
- Pros
 - Excellent temporal resolution
 - Allows visualization of arterial, capillary, and venous phases in a single acquisition
- Limitations
 - Lower spatial resolution than conventional single-phase CE-MRA
 - Requires gadolinium contrast





Sag TRICKS (mask)



SUB:Sag TRICKS (composite)



Subtraction coronal MPR



Early filling at T6 feeding cluster of vessels in the spinal canal, confirmed on DSA

UCSD Neuroradiology MRI Protocols

[Table of contents](#)

▼ MRI NEURO

Protocol	Brain	MR Angiogram	Head and Neck	Spine	Show All	
Brain MRI	Brain W/O	Brain +C	Brain Mets	Brain Tumor	Brain VWI +C	Ant Skull Base + Sella
	SRS/Brainlab	Pre-DBS	Spectroscopy			
	MS W/O	MS W/O LQ	MS +C	MS +C LQ		
	Epilepsy W/O	Epilepsy +C				
	Brain WO ARIA	Neuro Quant W/O	Neuro Quant +C			
	Neonate HIE +Spectroscopy	Neonate		Fetal		
	Neurocritical Care	Acute Stroke				
	CSF Flow - NPH	CSF Flow - Chiari	CSF Flow - Shunt	fMRI Guo/Hawley		
	Custom					
	MR Angiogram	Head MRA W/O	Head MRA +C	Custom MRA	VWI Head MRA W/O	
	Head MRV W/O	Head MRV +C				
	Neck MRA W/O	Neck MRA +C	Spine MRA			
Head and Neck MRI	Face W/O	FACE + C + Brain Lab	Neck W/O	Neck + C	Sella W/O	Sella + C
	Orbit W/O	Orbit + C	IAC W/O	IAC + C	Trigeminal W/O	Trigeminal + C
	Glossopharyngeal+C	Custom				
Spine MRI	C-Spine W/O	C-Spine + C	C-Spine Trauma	C-Spine CCJ	C-Spine MS W/O	
	C-Spine MS + C					
	T-Spine W/O	T-Spine +C	T-Spine Trauma	T-Spine MS W/O		
	T-Spine MS +C					
	L-Spine W/O	L-Spine +C	L-Spine Trauma	Spine Survey		

[Back to protocols](#)

Brain



▼ MRI NEURO

Protocol	Brain	MR Angiogram	Head and Neck	Spine	Show All	
Brain MRI	Brain W/O	Brain +C	Brain Mets	Brain Tumor	Brain VWI +C	Ant Skull Base + Sella
	SRS/Brainlab	Pre-DBS	Spectroscopy			
--	MS W/O	MS W/O LQ	MS +C	MS +C LQ		
--	Epilepsy W/O	Epilepsy +C				
--	Brain WO ARIA	Neuro Quant W/O	Neuro Quant +C			
--	Neonate HIE +Spectroscopy	Neonate	Fetal			
--	Neurocritical Care	Acute Stroke				
--	CSF Flow - NPH	CSF Flow - Chiari	CSF Flow - Shunt	fMRI Guo/Hawley		
--	Custom					

[Back to protocols](#)

Brain WO


- Indication:
 - New focal neurologic defect, fixed or worsening (suspected stroke)
 - Altered mental status
 - Dementia
 - Psychiatric disorder
 - Headaches
 - Symptomatic subacute or chronic head trauma
- Sequences:
 1. Sag T1 (FSE or FLAIR)
 2. Ax T1 (FSE or FLAIR)
 3. Ax T2 FSE
 4. Ax T2 FLAIR FS
 5. Cor T2 FSE
 6. Ax DWI / ADC / B0
 7. Ax 3D SWAN
 8. Ax 3D ASL




Acute Stroke

Stroke MRI Protocol: DWI, T2*GRE/SWAN, FLAIR and ASL perfusion imaging. MRA head or MRA Neck only if patient did not have CTA and is stable to continue. No additional sequences will be added. If the patient is moving during the DWI, the exam will not continue. *We do not do contrast MR perfusion for stroke codes.*


- Indication:
 - Abbreviated MRI appropriate for patients in the 6-24 hour or wake-up stroke window. The following steps must be taken:
 - Stroke team must put in "Acute Stroke Brain MRI" and MRA orders (if vessel imaging not performed by CTA)
 - Head CT and CXR screen for metal
 - Stroke Leader Calls MRI to discuss ETA and logistics
 - Radiologist calls MRI to give them a heads up and that metal was screened for
 - No discussion between stroke team and radiology is required
 - If the patient is moving during the DWI, the exam will not continue
- Sequences:
 - Standard
 - Ax DWI / ADC / B0
 - Ax FLAIR FSE
 - Ax SWAN
 - Ax ASL
 - If no CTA vessel imaging obtained
 - MRA is obtained




DWI




ADC




B0



T2 FLAIR




SWAN




ASL

Brain C+


- Indication
 - Infection
 - TIA
 - Parenchymal hemorrhage of uncertain etiology
 - Ataxia
 - Headache
 - Worsening mental status despite clinical management of the suspected underlying cause (intoxication, medication-related, hypoglycemia, sepsis)
- Sequences:
 - [Brain WO protocol](#) PLUS
 - SAG T2 FLAIR CUBE
 - Axial T1 +C (FSE or FLAIR)
 - Coronal T1 +C (FSE or FLAIR)
 - Axial T1 MPRAGE +C



Ax T1 +C




Cor T1 +C




T1 MP-RAGE +C

Brain Mets

- Indication
 - Metastatic brain lesions and meningiomas
 - Do not use for primary brain tumors (use [brain tumor](#) protocol)
 - If primary brain tumor and assessing leptomeningeal involvement, brain mets can be appropriate
- Sequences:
 - [Brain WO protocol](#) PLUS
 - Axial T1 +C (FSE or FLAIR)
 - Coronal T1 +C (FSE or FLAIR)
 - 3D Sag T1 MPRAGE +C
 - 3D Sag T2 FLAIR CUBE +C




Ax T1 +C



Cor T1 +C



T2 FLAIR CUBE +C




T1 MP-RAGE +C

Leptomeningeal metastases
occult on other sequences


[Back to protocols](#)

Brain Tumor


- Description:
 - Used for primary brain tumors, eg GBM, meningioma, medulloblastoma, etc
 - Do not use for metastatic brain disease
- Sequences:
 - [Brain +C protocol](#)
 - MINUS
 - Cor T1 FS C+
 - PLUS
 - Sag T1 MP-RAGE (plus reformats)
 - Sag T2 FLAIR FS CUBE (noncontrast)
 - Perfusion C+
 - Perfusion Field Map C+
 - Diffusion tensor imaging




Sag T1 +MP-RAGE



Sag T2 FLAIR



DSC Perfusion C+




Brain Operative Planning (epic button pending)


- For brain masses in the immediate post-operative period (24-72 hours), a full “Brain Tumor” protocol is not necessary. There is no prebuilt “epic” protocol button for this, but you can choose “Brain Mets” sequence and add any additional sequences requested by ordering (eg, T2 CUBE)
- Example:
 - “Reason for exam: Brain mass, include pre-post-contrast T1 FSPGR, T2 CUBE sagittally acquired with axial and coronal recons, sagittal FLAIR, Brainlab and DTI.”
- Sequences:
 - Can use “Brain mets” and add other requested sequences (eg, T2 CUBE)

Epilepsy (WO or W/WO)


- Indication:
 - Epilepsy
 - New seizure
- Sequences:
 - [Brain WO protocol](#) PLUS
 - Sagittal T1 MP-RAGE
 - Coronal oblique T2 FSE
 - Coronal oblique T2 FLAIR
 - Optional
 - Axial and coronal T1 +C
 - Sagittal T1 MP-RAGE +C



MP-RAGE



Cor oblq T2 FSE




Cor oblq T2 FLAIR


Coronal oblique sequences should be thin section and perpendicular to the long axis of the hippocampus.

MS (WO or W/WO, +/- lesionquant)


- Indication:
 - Multiple sclerosis or other demyelinating disease
 - Whether contrast is needed is variable, usually only if provider orders it, most often during initial work up and/or clinical progression/exacerbation
 - Follow up exams are most often WO contrast
- Sequences:
 - [Brain WO protocol](#) PLUS
 - Sag T1 FSPGR
 - Sag T2 FLAIR CUBE
 - Axial T1 +C (FSE or FLAIR)
 - Cor T1 FS +C FSE
- Notes
 - Does not include 3D post-contrast MP-RAGE




MP-RAGE



Sag T2 FLAIR CUBE




Ax T1 FSE C+




Cor T1 C+ FSE

Ant Skull Base Brain + Sella

- Indication
 - Use for all “MRI Brain + Sella” orders, all pre-op or immediate post-op studies for sella or anterior skull base masses
 - Operative planning and neuronavigation MRI for complex sella and third ventricle lesions
 - Includes large field for endonasal approach planning and volumetric for neuronavigation
 - [Sella/Pituitary protocol](#) (head and neck section below) is superior for initial diagnosis of sella/pituitary abnormalities, which takes less time to perform and has option of adding dynamic post-contrast
- Sequences:
 - Whole brain FOV
 - [Brain +C protocol](#) PLUS
 - Sag T1 MP-RAGE (instead of FSE)
 - Sag T2 CUBE (no Cor T2 FSE)
 - Small FOV (including sella and sphenoid sinus)
 - Axial T1 FSE
 - Coronal T2 sella to include sphenoid sinuses
 - Sag T1 FS + C sella to include sphenoid sinuses
 - Cor T1FS + C sella to include sphenoid sinuses




Sag T1 MP-RAGE
(replaces Sag T1 FSE)




Sag T2 CUBE
(replaces Cor T2 FSE)




Ax T1 FSE



Sag T1 FS C+







Cor T2 FSE



Cor T1 FS C+

Brain Vessel Wall Imaging


- Indication:
 - Vasculitis
- Sequences: MRI + MRA of the head
 - [Brain WO protocol](#)
 - MINUS
 - Ax T1
 - Cor T2
 - PLUS
 - 3D T2 CUBE
 - 3D MDSE Pre (3.0T) or 3D T1 Cube FS POST (1.5T)
 - 3D MDSE C+ (3.0T) or 3D T1 Cube FS POST (1.5T)
 - MRA
 - Ax 3D TOF



3D MDSE

3D MDSE C+

MRA TOF



3D MDSE

3D MDSE C+

MRA TOF

Circumferential vessel wall enhancement

[Back to protocols](#)

Head and Neck

▼ MRI NEURO

Protocol

Brain MR Angiogram **Head and Neck** Spine Show All

Head and Neck MRI

Face W/O FACE + C + Brain Lab Neck W/O Neck + C Sella W/O Sella + C
Orbit W/O Orbit + C IAC W/O IAC + C Trigeminal W/O Trigeminal + C
Glossopharyngeal+C Custom

Face +/- BrainLab

- Indication
 - Bell's Palsy
- Coverage
 - AP: Nose through brainstem
 - CC: frontal sinuses through entire chin.
 - Use 20 cm FOV for face unless very large patient. Do not include entire neck.
 - Brain lab to cover parotid, face, and skull base
- Sequences:
 - Ax T1 PROP
 - Cor T1 PROP
 - SAG T1 PROP
 - Ax STIR PROP
 - Cor STIR PROP
 - Ax MUSE DWI
 - +C Ax T1 FS FLEX (Water)
 - +C Cor T1 FLEX (Water)
 - Ax +C 3D T1 FS ISO (Brain Lab with same coverage of the face and skull base)



Ax T1 prop



Sag T1




Cor T1 prop




Cor STIR prop



Ax STIR PROP




Ax MUSE DWI



Water: +C Ax T1 FS FLEX



Ax +C 3D T1 FS ISO




Water: +c Cor T1 FLEX

[Back to protocols](#)

Neck


- Indication
 - Tumor, infection
- Sequences:
 - Cor T1 PROP
 - Sag T1 PROP
 - Ax DWI MUSE
 - Ax T1 PROP
 - Cor STIR PROP
 - Cor T1 +C
 - Ax T1 +C
- Dixon method:
 - Ax T2 (Water + In-Phase)
 - Ax T1 (Water + In-Phase)
 - Cor T1 (Water)
 - Sag T1 +C (Water)
 - Ax T1 +C (Water + In-Phase)
 - Cor T1 +C (Water)




Cor T1 PROP




Sag T1 PROP




Cor STIR PROP




Water: Ax T2 FSE Flex




InPhase: Ax T2 FSE Flex



Ax T1 PROP




Ax MUSE DWI




Water: Cor LAVA T1

Neck (continued)


- Indication
 - Tumor, infection
- Sequences:
 - Cor T1 PROP
 - Sag T1 PROP
 - Ax DWI MUSE
 - Ax T1 PROP
 - Cor STIR PROP
 - Cor T1 +C
 - Ax T1 +C
- Dixon method:
 - Ax T2 (Water + In-Phase)
 - Ax T1 (Water + In-Phase)
 - Cor T1 (Water)
 - Sag T1 +C (Water)
 - Ax T1 +C (Water + In-Phase)
 - Cor T1 +C (Water)




Water: Ax LAVA T1 FS



InPhase: Ax LAVA T1 FS




Water: Ax LAVA T1 FS C+




InPhase: Ax LAVA T1 FS C+




Water: Sag LAVA T1 FS C+



Ax T1 FS C+



Cor T1 FS C+




Water: Cor LAVA T1 FS c+

Sella/Pituitary


- Indication
 - Pituitary dysfunction, Sellar or suprasellar mass
- Sequences
 - Whole Brain FOV
 - Ax DWI
 - AX T2 FSE
 - AX T1 SE +C
 - Small FOV
 - Sag T1
 - Cor T1 SE
 - Cor T2 FSE
 - Cor T1 +C
 - Sag T1 +C
 - Comments
 - If post-operative, small FOV contrast sequences need to be done with **fat saturation**




Ax T2 FSE




Ax T2 FSE




Ax T1 FLAIR C+




Cor T2 FSE




Sag T1 FSE



Sag T1 F SE **FS** C+



Cor T1 FSE




Cor T1 FSE **FS** C+

Orbits – update*

- Indication
 - Orbit pathology, visual deficit
- Sequence
 - Sag T1
 - Ax T1 PROP
 - Cor T1 PROP
 - Ax DWI (whole brain)
 - Ax DWI MUSE
 - Cor STIR
 - T2 CUBE
 - Water: Cor T1 +C (FS)
 - Water: Ax T1 FLEX +C (FS)
- Optional:
 - Post Contrast Sag T1 Fat Sat along optic nerve (each side separate)
 - Add T2 FSE/ T2 FLAIR whole brain if visual field deficit/cranial nerve deficits
- FOV
 - FOV through orbits to include brainstem
 - Coronal perpendicular to optic nerves
 - Axial parallel to optic nerves




Sag T1 FSE




DWI




Cor T2 STIR




Water: Cor T1 +C (FS)



T2 CUBE




Ax T1 FS FSE




Water: Ax T1 FLEX +C

IAC


- Indication:
 - Tinnitus, vertigo, hearing loss, CN6 or 7 palsy, 7th nerve palsy
 - Do NOT use for known schwannoma, sudden onset hearing loss, or presurgical evaluation (cochlear implant)
- Sequences
 - Whole Brain:
 - Ax DWI
 - Ax T1 FLAIR +C whole brain (no Fat Sat)
 - IAC
 - Sag T1
 - Ax DWI FSE PROP (non-echo planar DWI)
 - Cor T2
 - Ax T1
 - 3D T2 Cube or FIESTA
 - Ax T1 FS PROP +C
 - Cor T1 FS PROP +C
- Comments:
 - For cranial nerves extend slices anteriorly to posterior orbit & superiorly to 3rd ventricle




Sag T1 FSE




DWI




DWI PROP




Ax T1 FSE C+




Cor T2 FSE




Cor T1 FS FSE C+



Ax FIESTA




Ax T1 FSE




Ax T1 FS FSE C+

Trigeminal


- Indication:
 - Trigeminal neuralgia
- Sequences
 - Ax 3D TOF MRA
 - Ax 3D T2 Cube or FIESTA
 - Ax 3D MP-RAGE +C




Ax 3D SPGR



Ax 3D T2 Cube




Ax 3D SPGR C+



3D reconstructions


Glossopharyngeal

- Indication:
 - Glossopharyngeal Neuralgia
- Coverage
 - Smaller FOV: coverage top of pons through C3, brain stem through all of face/mouth
- Sequences
 - Whole brain
 - Ax DWI
 - Ax T2 FLAIR
 - Ax T1 FS +C
 - Small FOV
 - Sag T1
 - Time of flight MRA (does not include full circle of Willis)
 - 3D T2 Cube or FIESTA Cor post contrast fat saturated
 - Cor T1 FS +
 - T1 MP-RAGE C+




Head MRA

- Non-contrast is sufficient for intracranial aneurysm screening and LVO detection (eg, time of flight and phase-contrast techniques)
- Contrast useful for more accurate vessel lumen assessment (especially in the neck) for high-grade carotid stenosis, dissections, and aortic arch pathology, and can also be useful to characterize treated/coiled aneurysms to look for residual filling
- Sequences:
 - Ax 3D TOF FSPGR
 - SAG 3D Inhance Pre
 - +/- Ax 3D TOF C+ FSPGR




Ax 3D TOF



Ax 3D TOF C+


Head MRV

- For venous sinus thrombosis, contrast-enhanced is **strongly** preferred (push for this).
- Pregnant patients cannot receive contrast
- Sequences:
 - SAG 3D Inhance Pre
 - Cor 2DTOF SPGR
 - Sag 2DTOF SPGR (slight oblique angle)
 - Cor 3DTOF FSPGR +C




3D Sag Inhance Pre

3D Baseline



Cor 2D TOF




Sag 2D TOF

3D C+


[Back to protocols](#)

Neck MRA


- Non-contrast is sufficient for intracranial aneurysm screening and vessel occlusion (eg, time of flight and phase-contrast techniques)
- Contrast useful for more accurate vessel lumen assessment (especially in the neck) for high-grade carotid stenosis, dissections, and aortic arch pathology, and can also be useful to characterize treated/coiled aneurysms to look for residual filling
- FOV:
 - Aortic arch to circle of Willis
- Sequences:
 - Standard
 - InHance
 - Ax 3D TOF SPGR
 - Cor 3D TOF FSPGR+C
 - For Dissection/Trauma:
 - AX T1 FS GATED BLACK BLOOD through skull base



Cor 3D InHance (subtracted vs background)




Ax 3D TOF



Black Blood T1 FS or MDSE

Spine MRA


- Indication:
 - Spinal arteriovenous shunt lesions
- Sequences:
 - Sag TRICKS +C (mask)
 - Subtraction images/phases



Sag TRICKS (mask)



SUB:Sag TRICKS (composite)



Subtraction coronal MPR

Spine

MRI NEURO	
Protocol	Brain MR Angiogram Head and Neck Spine Show All
Spine MRI	C-Spine W/O C-Spine + C C-Spine Trauma C-Spine CCJ C-Spine MS W/O
	C-Spine MS + C
--	T-Spine W/O T-Spine +C T-Spine Trauma T-Spine MS W/O
	T-Spine MS +C
--	L-Spine W/O L-Spine +C L-Spine Trauma Spine Survey

C-Spine WO

- Indication:
 - Disc disease, degenerative, pain, radiculopathy, including post-op
- Sequences:
 - Standard
 - Sag T2 FSE
 - Sag T1 FSE
 - Sag STIR
 - Ax T2 FSE
 - Ax MERGE
 - Optional
 - Cor T1 FSE
 - Cor T2 FSE
- Comments
 - Axial imaging from C2 to top of T1
 - For scoliosis and tethered cord add coronal




Sag T1 TSE




Sag T2 FSE



Sag T2 STIR




Ax T2 FSE




3D Ax <MERGE

C-Spine CCJ (craniocervical junction)


- Indication:
 - Trauma, craniocervical ligament evaluation
- Coverage:
 - Axials go to sphenoid sinus to T1-T2 (higher than standard C-spine WO)
- Sequences:
 - C-Spine Trauma WO PLUS
 - Sag 3D T2 Cube
 - Cor STIR
 - Cor PD
 - Small FOV centered at C1-2
 - 3D PD Cube
 - Cor STIR



Sag T2 CUBE



Ax T2 FSE FS







Cor T2 STIR



3D PD cube


C-Spine CCJ (craniocervical junction)

Partial tear of the left **alar ligament**
and **apical ligament**



C-Spine W

- Notes:
 - Indications: Tumor, Infection, Syring
- Sequences:
 - For post-contrast, some scanners utilize FLEX (Dixon), others FSE
 - [C-Spine WO PLUS](#)
 - Sag T1 FS + C
 - Axial T1 FS +C
 - If post fat saturation is suboptimal, repeat sequences non fs.
- Optional
 - For scoliosis, tethered cord and Neurofibromatosis, add coronal
 - If there is concern for syrinx, add 1-mm 3D sagittal FIESTA through the spinal cord.




Water: Ax T1 FS FLEX +C




Water: Sag T1 FS FLEX +C

C-Spine MS (WO or WO/W)

- All scanners utilize FSE (not FLEX)
 - T1 and T2 FSE sequences maintain more consistent T1 and T2 contrast ensuring optimal cord-to-lesion contrast, and because they are not influenced by FLEX artifacts (Dixon fat-water separation processing), which could mimic or obscure MS plaques.
 - Maintains consistency with prior studies, which is important for accurately assessing new vs. chronic lesions.
- Sequences:
 - [C-Spine WO PLUS](#)
 - Sag PD
 - Sag T1 +C (no Fat Sat)
 - Axial T1 +C (no Fat Sat)



Sag T1 PD




Ax NOFS T1 C+

Sag NOFS T1 C+

[Back to protocols](#)

T-Spine WO

- Indication:
 - Disc disease, pain, radiculopathy
- FOV
 - Axial coverage through ROI
 - Angle sagittals to spine, if scoliosis, do upper and lower angles.
- Sequences:
 - Standard
 - Sag T1 FSE counter
 - Sag T1 FSE or FLAIR
 - Sag T2 FSE
 - Sag STIR
 - Ax MERGE
 - Ax T2 FSE
 - Optional
 - Cor T1 FSE
 - Cor T2 FSE
 - For scoliosis and tethered cord, add coronal



Sag T1 counter




Sag T1 FLAIR




Sag T2




Sag STIR




Ax 2D MERGE upper



Ax T2 FSE upper



Ax 2D MERGE lower



Ax T2 FSE lower

T-Spine Trauma

- Indication:
 - Trauma
 - Acute thoracic or lumbar spine injury detected on radiographs or noncontrast CT. Neurologic abnormalities.
- Sequences:
 - [T-Spine WO](#)
 - Sag T1 FSE
 - Sag T2 FSE
 - Sag STIR
 - Sag GRE
 - Sag T2 CUBE
 - Ax T2 FSE FS, upper and lower
 - Optional
 - Ax GRE
 - Cor T1 FSE
 - Axial coverage through ROI
 - For scoliosis and tethered cord add coronal
 - Angle sagittals to spine, if scoliosis, do upper and lower angles.



Sag counter




Sag T2 FSE




Sag T1 FSE




Sag STIR




Sag GRE



Sag T2 CUBE





Ax T2 FSE FS upper



Ax T2 FSE FS lower

T-Spine WO/W

- Indication:
 - Tumor, Infection, Syrinx
- Sequences:
 - [T-Spine WO](#) PLUS
 - Sag T1 +C FSE FS
 - Ax T1 +C FSE FS
 - If Fat Sat fails, do non Fat Sat post contrast
- Optional
 - Cor T1 FSE
 - Cor T2 FSE
 - For neurofibromatosis, add coronal T2



T-Spine MS

- Indication:
 - MS
- Sequences:
 - Sag STIR
 - Sag PD
 - Sag T2 FSE
 - Sag T1 FS
 - Ax MERGE
 - Ax T2 FSE
 - Sag T1 +C (no Fat Sat)
 - Axial T1 +C (no Fat Sat)
- Angle sagittals to spine, if scoliosis, do upper and lower angles.



Sag T2 count



Sag T2



Sag T1




Sag PD




Sag STIR




Sag T1 No FS C+




Ax T2 FSE upper





Ax MERGE upper



Ax T2 FSE lower





Ax MERGE lower



Ax T1 upper and lower


L-Spine WO

- Indication:
 - Disc disease, pain, radiculopathy
- Sequences:
 - Standard
 - Sag T1 FSE
 - Sag T2 FSE
 - Sag STIR
 - Ax T2 FSE
 - Ax T1 FSE
 - Optional (eg, scoliosis, tethered cord)
 - Cor T1 FSE
 - Cor T2 FSE





L-Spine Trauma

- Indication:
 - Trauma
 - Acute thoracic or lumbar spine injury detected on radiographs or noncontrast CT. Neurologic abnormalities.
- Sequences:
 - [Lumbar Spine WO](#) PLUS
 - Sag GRE
 - Ax T2 FSE FS (in place of non-FS)
 - Optional
 - Cor T1 FSE
 - Cor T2 FSE



L-Spine C+

- Indication:
 - Tumor, post-operative, infection
- Sequences:
 - [L-Spine WO](#) PLUS:
 - Sag T1 +C FSE FS
 - Axial T1 FS +C
 - Optional
 - Cor T1 FSE
 - Cor T2 FSE
 - If metal hardware in spine, obtain the post contrast without fat sat
 - For neurofibromatosis add coronal T2



Spine Survey

- Indication:
 - DO NOT USE FOR MS OR TRANSVERSE MYELITIS
 - Metastases
 - Non-localized infection
 - Acute myelopathy / cord compression
 - [See next slide for additional details](#)
- FOV
 - 2 sets of sagittal images, upper (cervicothoracic) and lower (thoracolumbar)
- Sequences:
 - Water: Sag T2 FLEX
 - InPhase: Sag T2 FLEX
 - Sag T1 PROP
 - Water: +C Sag T1 FLEX
 - InPhase: +C Sag T1 FLEX
 - C+ Ax 3D LAVA (upper, thoracolumbar, lower)
- Comments
 - For scoliosis, tethered cord or neurofibromatosis, add coronal
 - Post-contrast images in patient with spinal hardware need to be done non-fat suppressed, and have axial T1 pre-contrast for comparison




Water: Sag T2 FLEX



InPhase: Sag T2 FLEX




Sag T1 PROP




Water: +C Sag T1 FLEX



InPhase: +C Sag T1 FLEX



Water: Sag T2 FLEX



InPhase: Sag T2 FLEX




Sag T1 PROP




Water: +C Sag T1 FLEX




InPhase: +C Sag T1 FLEX



C+ Ax 3D LAVA Upper



C+ Ax 3D LAVA Thoracolumbar



C+ Ax 3D LAVA Lower

[Back to protocols](#)




Spine Survey

- Advantages:
 - Faster
 - Patients tolerate the exam
 - Not tailored for the evaluation of degenerative changes
- Disadvantages:
 - Larger FOV decreases lesions conspicuity
 - Lack of axial T2 – limited evaluation of nerve roots
 - Not tailored for the evaluation of degenerative changes
 - Surgical or radiation planning may require dedicated exam of spine segment

Water: +C Sag T1 FLEX




C+ Ax 3D LAVA



Cases

Case 1



Order: Neck MRA

Reason for exam: Dissection

What sequence should we add?

- A. Time of Flight
- B. Black blood fat-saturated
- C. T1 non-fat saturated MRA
- D. Phase contrast

Order: Neck MRA





Reason for exam: Dissection

What sequence should we add?

- A. Time of Flight
- B. Black blood fat-saturated
- C. T1 non-fat saturated MRA
- D. Phase contrast

B. Black blood fat-saturated

Black blood: T1 Fat-saturated



Reason for exam: Aneurysm

- A. Brain C+
- B. MRA Head WO
- C. MRA Head C+
- D. Change to CTA

Reason for exam: Aneurysm

- A. Brain C+
- B. MRA Head WO
- C. MRA Head C+
- D. Change to CTA

B. MRA Head WO

Reason for exam: Aneurysm s/p coiling

A. MRA


B. CTA

Reason for exam: Aneurysm s/p coiling

A. MRA

B. CTA

A. MRA



Reason for exam: Aneurysm s/p clipping

- A. MRA Head C+
- B. CTA
- C. MRA Head WO


Reason for exam: Aneurysm s/p clipping

A. MRA Head C+

B. CTA

C. MRA Head WO

C. CTA






Aneurysm coiling

- MRA is superior to CTA
- Significant beam-hardening and streak artifact on CTA
- Contrast can increase visualization of residual aneurysm
- Exception: Flow-diverter stent, CTA is more useful for evaluation in-stent lumen/stenosis

Aneurysm clipping


- CTA is superior to MRA
- Clips cause marked metallic susceptibility artifact at MRI
- CTA is first-line imaging modality



MR Angiography


- Most techniques use T1-weighted echo gradient sequences
- Most suppress background signal of stationary tissue (except black blood)
- All techniques can be adapted to 3D

Time of Flight (non-contrast)




Adapted from <https://doi.org/10.1148/rg.2021210141>

Saturation-band for incoming venous signal




2D TOF

- Fine slices are stacked alongside each other
- Sensitive to slow flows (which do not remain in the slice for long and will therefore not be saturated)
- Poor spatial resolution along the axis of the slice stack




3d TOF


- Good spatial resolution in all 3 planes
- Better signal-to-noise ratio
- Longer acquisition



2D TOF



3D TOF



Phase-Contrast MRA

Signal in arteries + veins

Axial



Phase-Contrast (no gadolinium)

Bipolar gradient
(VENC)


Stationary Spin




Flowing Spin




Challenge case



3D InHance
(phase contrast)





3D TOF



CTA


Subclavian Steal




Challenge case



3D InHance
(phase contrast)






2D TOF



MRA C+

Jugular venous reflux





Reason for exam: Cerebral venous thrombosis

- A. MRA C+
- B. MRV WO
- C. MRV C+
- D. MRI Brain WO/W

Reason for exam: Cerebral venous thrombosis

- A. MRA C+
- B. MRV WO
- C. MRV C+
- D. MRI Brain WO/W

C. MRV C+



Reason for exam: Acute focal neurological deficit, stroke


- A. Acute stroke
- B. Brain WO
- C. Brain C+
- D. MRA Head

Reason for exam: Acute focal neurological deficit, stroke


- A. Acute stroke
- B. Brain WO
- C. Brain C+
- D. MRA Head

B. Brain WO


Brain WO




T1




T1




T2




T2 FLAIR




T2




DWI
(B 1000)




B0




ADC



SWAN



Phase map




ASL


Acute Stroke

e HIE + Spectroscopy	Neonate
ritical Care	Acute Stroke
ow - NPH	CSE Flow - Chiari
	CSE Flo


- For patients who are candidates for intervention (in the 6-24 hour or wake-up stroke window)
- Not for “standard” stroke work-up
- **Stroke team** must put in “Acute Stroke Brain MRI” orders




DWI




ADC




B0



T2 FLAIR




SWAN




ASL

Which sequence in the routine Brain WO protocol is the fastest sequence to perform?





DWI
~1 min



Diffusion weighted imaging


Single-shot echo planar imaging (standard)

- Acquires all of k-space lines with one RF excitation pulse (single TR) per slice
- Fast acquisition (few motion artifacts)
- Prone to susceptibility artifacts (distortion near air-bone interfaces)



MUSE (Multiplexed Sensitivity-Encoding)

- Multi-shot EPI sequence designed to improve resolution and reduce distortion compared to standard single-shot EP
- Acquires data over multiple shots rather than a single shot
- Higher resolution
- Reduced geometric distortion (better for regions affected by susceptibility artifacts, such as the skull base and orbits)
- Longer scan times, more motion artifact



DWI MUSE


DWI PROP



PROPELLER


- Non-EPI DWI (utilizes fast-spin echo), acquires multiple overlapping radial "blades" in k-space and averages them for artifact reduction.
- Significantly reduced motion artifacts
- Lower susceptibility artifacts than ss-EPI (less distortion in air-bone interfaces like the skull base)
- Higher SNR than standard DWI due to multiple overlapping k-space acquisitions.
- Longer scan times

Diffusion weighted imaging



B-1000


÷



B0




ADC



Quantitative



Even more math



eADC

Qualitative

ADC tells you whether diffusion is truly restricted;
eADC just makes restricted areas easier to see

Reason for exam: Stroke, “thin cuts” through the brain stem.

What additional sequence would probably be most helpful for posterior circulation infarct?


- A. Coronal DWI
- B. Axial thin slice DWI
- C. T2 CUBE IAC
- D. FIESTA IAC

Reason for exam: Stroke, “thin cuts” through the brain stem.


What additional sequence would probably be most helpful for posterior circulation infarct?

- A. Coronal DWI
- B. Axial thin slice DWI
- C. T2 CUBE IAC
- D. FIESTA IAC

A. Coronal DWI




DWI in a different plane (eg, coronal) is often most helpful (if had to pick one)



FIESTA

- Slightly higher **spatial resolution**
- More CSF pulsation artifact
- Faster



FIESTA



T2 CUBE

- Better T2 weighting and **tissue contrast**
- Less CSF pulsation artifact
- Slower



T2 CUBE



Reason for exam: Headaches, blurred vision, evaluate for pseudotumor cerebri

- A. Brain WO
- B. Brain C+
- C. Brain Mets
- D. Head MRV WO

Brain C+



Which of the following pathologies should NOT receive a “Brain Tumor” protocol?

- A. Meningioma
- B. Glioblastoma
- C. Medulloblastoma
- D. Hemangioblastoma

Which of the following pathologies should NOT receive a “Brain Tumor” protocol?

- A. Meningioma
- B. Glioblastoma
- C. Medulloblastoma
- D. Hemangioblastoma

A. Meningioma

Reason for exam: Meningioma

- A. Brain Mets
- B. Brain Tumor
- C. Brain C+
- D. Head MRA

Reason for exam: Meningioma

- A. Brain Mets
- B. Brain Tumor
- C. Brain C+
- D. Head MRA

C. Brain C+

Brain C+

✓ MRI NEURO

Protocol

Brain

MR Angiogram

Head and Neck

Spine

Show All

Brain MRI

Brain W/O

Brain +C


Brain Mets

Brain Tumor


Brain VWI +C

Ant Skull Base + Sella


- Brain WO protocol PLUS:



Cor T1 FS +C



Ax T1 +C



T1 MP-RAGE +C

Brain Mets

✓ MRI NEURO

Protocol

Brain

MR Angiogram

Head and Neck

Spine

Show All

Brain MRI

Brain W/O

Brain +C

Brain Mets

Brain Tumor


Brain VWI +C

Ant Skull Base + Sella


- Brain WO PLUS:




T2 FLAIR CUBE +C



Cor T1 FS +C



Ax T1 +C



T1 MP-RAGE +C

T2 FLAIR shows contrast enhancement due to mild T1 weighting properties from its long inversion recovery pulse, despite being heavily T2-weighted

[Table of contents](#)

Brain Tumor

✓ MRI NEURO

Protocol

Brain

MR Angiogram

Head and Neck

Spine

Show All

Brain MRI

Brain W/O

Brain +C

Brain Mets

Brain Tumor


Brain VWI +C

Ant Skull Base + Sella


- Brain WO PLUS:




T2 FLAIR CUBE +C




Cor T1 FS +C



Ax T1 +C




T1 MP-RAGE +C




DSC Perfusion C+

Post-contrast Brain protocols


Brain C+ Default



Cor T1 FS +C




Ax T1 +C




T1 MP-RAGE +C


Brain Mets Leptomeningeal
(infection, mets)




Cor T1 FS +C



Ax T1 +C




T1 MP-RAGE +C




T2 FLAIR CUBE +C


Brain Tumor Primary* brain tumor
(glioma)



Ax T1 +C



T1 MP-RAGE +C



DSC Perfusion C+

[Table of contents](#)

Which of the following is NOT an indication for MRI spine survey?

- A. Prostate Cancer
- B. Trauma
- C. Non-localized infection
- D. Drop metastases


Which of the following is NOT an indication for MRI spine survey?

- A. Prostate Cancer
- B. Trauma
- C. Non-localized infection
- D. Drop metastases


B. Trauma



Water: Sag T2 FLEX




Water: +C Sag T1 FLEX



Water: Sag T2 FLEX



Water: +C Sag T1 FLEX



C+ Ax 3D LAVA Thoracolumbar

Spine survey has large FOV and does not provide axial T2/T1 precontrast
Contrast unnecessary for trauma

Reason for exam: AMS with first time seizure

- A. Brain WO
- B. Brain WO/W
- C. Epilepsy C+
- D. More info

Reason for exam: AMS with first time seizure


- A. Brain WO
- B. Brain WO/W
- C. Epilepsy C+
- D. More info

B. Brain WO/W

Epilepsy imaging is best at high magnet strength 3T – no inpatient 3T scanner at JMC

Epilepsy imaging requires high-quality images and patient remaining still

Epilepsy protocol does NOT come with MP-RAGE






Reason for exam: Cervical spine fracture, C2

- A. Cervical spine WO
- B. Cervical spine CCJ
- C. Spine Survey
- D. Cervical Spine MS

Reason for exam: Cervical spine fracture, C2

- A. Cervical spine WO
- B. Cervical spine CCJ
- C. Spine Survey
- D. Cervical Spine MS

B. Cervical Spine CCJ



Reason for exam: Status post EEA, Dr. Beaumont protocol

- A. Ant Skull Base + Sella
- B. Brain C+
- C. Maxillofacial C+
- D. IAC C+

Reason for exam: Status post EEA, Dr. Beaumont protocol


- A. Ant Skull Base + Sella
- B. Brain C+
- C. Maxillofacial C+
- D. IAC C+

A. Ant Skull Base + Sella


Reason for exam: Lesion localizing to cavernous sinus

- A. CT venogram
- B. MRI Sella WO/W
- C. MRI Orbits WO/W
- D. MRI Brain Mets


C. MRI Orbits WO/W (OK to add MRI brain mets as it includes ASL)




DWI MUSE



T1 C+



T2



Lesion nearly invisible on MRI Brain MP-RAGE

Reason for exam: MRA with multifocal intracranial arterial narrowing.
Concern for intracranial vasculitis.


- A. CTA Head and Neck
- B. MRI Brain C+
- C. MR Vessel Wall Imaging (VWI)
- D. MR Perfusion

Reason for exam: MRA with multifocal intracranial arterial narrowing.
Concern for intracranial vasculitis.


- A. CTA Head and Neck
- B. MRI Brain C+
- C. MR Vessel Wall Imaging (VWI)
- D. MR Perfusion

C. MR Vessel Wall Imaging (VWI)

Vessel Wall Imaging Protocol



3D MDSE




3D MDSE C+




MRA TOF

Enhancement patterns:

- Vasculitis – Concentric with positive remodeling
- ICAD – Eccentric
- Moyamoya – Concentric with negative remodeling
- RCVS – No enhancement



3D MDSE



3D MDSE C+



MRA TOF